



5358, 5345, 5394



5338-100



5338-70, 5393

## SPECIFICATIONS

### Headrest Translation

- Anterior Range ..... 5 in.
- Longitudinal Range ..... 8.5 in.
- Lateral Tilt ..... 30°
- Patient Weight Capacity ..... 400 lbs.
- Ultra™ Weight ..... 29 lbs.

## ORDERING INFORMATION

- 5370 — Ultra Shoulder Positioner™ and Cart.
- 5358 — Ultra Shoulder Positioner™  
Comes with two Patient Care Kits
- 5345 — Patient Care Kit  
Three per case (Others available)
- 5394 — Aluminum Side Rail Socket  
(Two required)
- 5338-24 — Shoulder Positioner Cart  
This specially designed cart transports the  
positioner, accessories and reserve Patient  
Care Kits
- 5338-100 — Lateral Support Pad Assembly  
This specialty pad assembly provides  
additional support and prevents the patient  
from shifting laterally during the rigors of  
surgery (Clark Socket included)
- 5338-70 — Translating Arm Board for may be  
positioned to provide a platform for either arm
- 5393 — Clark Socket

Attaches to any OR table.

**ULTRA SHOULDER POSITIONER™**  
FOR ARTHROSCOPIC AND OPEN SHOULDER SURGERY



MIZUHO OSI  
30031 AHERN AVENUE  
UNION CITY, CA 94587  
(510) 429-1500  
(800) 777-4674  
FAX (888) 329-4674  
FAX (510) 429-8500

MIZUHOSI.COM  
NEWHIPNEWS.COM  
NEWSPINENEWS.COM

ISO certified by TUV Rheinland of North America



Note: Mizuho OSI is constantly improving its products, all specifications are subject to change without notice.

Mizuho OSI is a Delaware Corporation. Manufactured in the USA.

©Mizuho OSI NW0377 Rev. C. 11/09



# ULTRA SHOULDER POSITIONER™

The Ultra Shoulder Positioner™ provides a safe, convenient and stable method of maintaining patients in the “beach chair” position for arthroscopic open shoulder surgery and joint replacement.

- Minimizes setup tasks
- Eliminates lateral decubitus positioning concerns
- Fits all surgical table models
- Offers effective solution for precise shoulder positioning



Headrest wings retract posteriorly creating a totally obstruction-free space around the patient's head for the anesthetist.



With the wings in place and head strap positioned, the headrest becomes a secure head support.

The cervical flexion angle for a kyphotic patient is easily accommodated with the adjustable headrest.



## ENHANCED EASE OF USE

The Ultra is designed to suit the high demands of the OR Suite.

- The center mounted support rod assures unobstructed radiolucency at the surgical site.
- The radiolucent backrest has a low contour for optimal shoulder access both anteriorly and posteriorly.
- All release levers and adjustment knobs are readily located for convenient access and easier positioning of the patient.
- The headrest adjusts superiorly and inferiorly to accommodate the patient's height.

## PATIENT CARE KITS

Single use Patient Care Kits may be customized to accommodate each patient. On the operative side, tear along the perforations to increase shoulder access. Leave pad in place to support the non-operative shoulder.



The headrest tilts laterally up to 30° in either direction.